

Torvan Medical Endosink

Purpose built to meet the needs of endoscope reprocessing workflows our Endosink provides highly customizable, ergonomic features within a durable station.



Overhead Tilt Adjust LED Lighting



Pegboard Storage & Accessory Area



Wrist Friendly Edges & Ergonomic Features



Programmable Height Adjustment



Bottom Filling Sink Basins with Etched Gallon Measures & Ergonomic External Pull Drains



BUILT FOR IMPROVED PATIENT CARE IN BUSY MEDICAL ENVIRONMENTS

Built for heavy daily use in multi-user environments the programmable height adjustment, overhead tilt-adjust lighting and vertical pegboard storage area allows versatile and ergonomically optimized placement of all storage and accessories within easy reach for all users.

PRODUCT	PRODUCT NUMBERS	FOOTPRINT DIMENSIONS	SURFACE HEIGHT	INCLUDED FEATURES
Endosink	T2-ENDOSINK-HA	48" W x 30" D	34"-45" Programmable Adjustable Height <small>*Fixed Height Available on Request</small>	<ul style="list-style-type: none"> • Two 20" x 20" x 8" Basins w. Gallon Etching • LED Overhead Tilt Light • Ergonomic Externally Operated Pull Drains • Backsplash Mount 10" Swing Spout Faucet • 2x Bottom Fill Mixers & Basin Outlets • Lock & Leveling Swivel Casters • Marine & Wrist Friendly Edges • 24" High Modular Pegboard Accessory Area

*Options and Features May Change



PROUDLY BUILT IN CANADA
TORVANMEDICAL.COM

